

This Page Is Inserted by IFW Operations
and is not a part of the Official Record

BEST AVAILABLE IMAGES

Defective images within this document are accurate representations of the original documents submitted by the applicant.

Defects in the images may include (but are not limited to):

- BLACK BORDERS
- TEXT CUT OFF AT TOP, BOTTOM OR SIDES
- FADED TEXT
- ILLEGIBLE TEXT
- SKEWED/SLANTED IMAGES
- COLORED PHOTOS
- BLACK OR VERY BLACK AND WHITE DARK PHOTOS
- GRAY SCALE DOCUMENTS

IMAGES ARE BEST AVAILABLE COPY.

**As rescanning documents *will not* correct images,
please do not report the images to the
Image Problems Mailbox.**

09735519.121300




Fig 1

Study Setup

- Clinical Study Definition
- Describes Sponsors & Investigators
- Declares Subject Attributes to Capture
- Associates Specific Lab Procedures with a Clinical Study
- Defines Genotype Results to Report

Study Protocol		Sponsor		Investigator	
Study Number 632-001	Sponsor Study Number 199-005	Sponsor Name A CO, JOHN DOE, Ph.D.	Investigator Name JOHN DOE, Ph.D.	Sponsor Company A CO.	Investigator Company A CO.
Sample Test PS1	Sample Test SS1	Sponsor Representative DAVID JONES	Investigator Representative JOHN DOE, Ph.D.	Sponsor Representative DAVID JONES	Investigator Representative JOHN DOE, Ph.D.
Sponsor Protocol Title Pharmacokinetic Evaluation of Olanzapine (olanzapine chloride) and IR cypofluorin Administered Alone and in the Presence of Ketonecarboxide		Investigator Protocol Title Genomic DNA isolation and molecular genotyping analysis of CYP2D6			
Procedures [HIA] Isolation, 3 ml. whole blood, Purpure K2 SpectroMax DNA quantitation CYP2C93 CYP2C92 Ver. 7		Subject Attributes Subject Number Gender Birthdate Ethnicity			
Created 10/4/99 12:02		Modified 10/5/99 14:38		Define Results	
DUALS		DUALS		New Modify Delete Save Cancel Done	

FIG. 2

Clinical Sample

Registration

- Provides Validation Checks for Accession & Tube ID's
- Accommodates Multiple Sample Tubes
- Enforces Controlled Subject Attribute Terms
- Supports Sample Workflow

FIG. 3

FIG. 3

Sample Tracking

- Supports Multiple Container Classes
- Allows User Defined Container Geometries & Templates
- Maintains Sample & Container Location
- Permits Flexible Sample Loading & Rearrangement
- Tracks and Maintains Container & Sample Ownership




FIG. 4

- Automates Laboratory Process Tracking
- Supports Standard Operating Procedures
- Maintains Uniform Laboratory Processes
- Records Chain of Custody
- Tracks Repeat Operations

FIG. 6

FIG. 6

Procedure		Status	SOP Number	SOP Version
DNA Isolation, 3 mL whole blood, Pargene Kit		APPROVED	GEN9709	C
2D6 Allele 'A' Identification		APPROVED	CYP2D6A	A
SpectroMax DNA quantitation		APPROVED	MAX9802	A
CYP2C9-3		APPROVED	CYP2C9-3	A
CYP2C9-2 Ver. 7		APPROVED	CYP2C9-2	A
CYP2C9-2 Ver. 6		APPROVED	CYP2C9-2	A
Procedure Name				
2D6 Allele 'A' Identification				
Procedure Description				
Laboratory Protocol for Identification of CYP2D6 'A' Allele by AFLP-PCR				
Gene		SOP Number		
[M]		CYP2D6A		
Alleles		SOP Version		
[Image]		A		
Status		APPROVED		
New		Modify		
Delete		Close		
Created		Modified		
DNALIMS		DNALIMS		
10/4/99 12:02		10/5/99 10:59		
Print		Save As		

FIG. 7

Procedure Steps

- A Single Step in a Lab Procedure
- Multiple Types:
 - Transfer
 - Dilution
 - Concentration Adjustment
 - Sample Preparation
- Highly Customizable
- Plug-in Architecture to Add New Types
- Interfaces to Automation

Procedure Steps

Procedures: DNA Isolation, 3 mL whole blood, Purgene Kit

Step	Step Input	Type	Functional Type
Thaw frozen blood		CheckBox	Batch
Gently mix sample		CheckBox	Batch
Transfer 3 mL of blood to Lysis tube		Functional	Batch
Add 9 mL of RBC lysis to RBC lysis tube		Functional	Batch
Mix and incubate 10 minutes at room temperature		CheckBox	Batch
Centrifuge 10 minutes at 3000 RPM		Text	Batch
Pour off supernatant into bechardous waste cont.		CheckBox	Batch
Resuspend cell pellet by vortexing		CheckBox	Batch
Add 3 mL of Cell Lysis Solution		Functional	Batch
Sample can be stored for 18 months at RT in Cell		Informational	Batch
STMP POINT		Informational	Batch

Step Type

☒ Informational ☒ Functional ☒ Sample

Reagent Addition

Volume: 9000.2371
 To final volume: FALSE
 Reagent name: RBC Lysis S
 Reagent prefix: RBC
 Volume optional: TRUE
 Wave scanning: TRUE
 Lock parameter: TRUE

Save Sequence

New Save
 Delete Cancel
 Modify Close

Created: 10/5/99 09:00 Modified: 9/30/99 14:42
 DIALIMS

FIG. 8

Genotype Results

- Accommodates Values for Multiple Genes, Alleles & Assays
- Provides Master Review by Accession Number
- Supports Acceptance & Final Approval by Study Director
- Imports Results Electronically or Manually

Study Protocol		Accession No.		Status
PS1	PS1	A3		OPEN
PS1	PS1	A4		OPEN
PS1	PS1	A5		OPEN
PS1	PS1	A6		OPEN
PS1	PS1	A7		OPEN
PS1	PS1	A8		OPEN
PS1	PS1	A9		OPEN

Accession No.		Status
AG		OPEN

Genotyping		DNA Purification		Chain Of Custody	
GENE	STATUS	INTERPRETATION	EXTENSIVE METABOLIZER	ENTERED BY	ENTERED ON
CYP2C9	COMPLETE			DNALIMS	10/6/99 6:48

Gene		Batch		Procedure		Created By		Created On	
m2	Alb	m/m		CYP2C9 Ver. 7	DNALIMS	10/6/99 6:09 48 PM	DNALIMS	10/6/99 6:09 48 PM	DNALIMS
m2	Alb	m/m		CYP2C9 Ver. 7	DNALIMS	10/6/99 6:09 48 PM	DNALIMS	10/6/99 6:09 48 PM	DNALIMS
FINAL RESULT									

FIG. 9

- Study
- Lab Procedures
- Sample
- Results

Flexible Audit Reporting Chain of Custody by Accession

FIG. 10

Reporting

Multiple Report Types

- Genotype Results
- DNA Integrity
- Purification Results
- Sample Lists
- Audit Trails

Flexible Reporting Output Using Excel

Customized Reporting Using 3rd Party Tools

DNA Purification Results										Sponsor:	
Study:										Representative:	
Test										Investigator(s):	
Study #: 999-xxx											
Sample #	Storage #	µg/ml	Protocol	Start Vol.	DNA Vol.	Yield	A260	A230	A280	Subject Number	
EA111111	0	Sample test	09/06/00	499	0	0	0	0	0	Initial date of Birth	
EA111112	0	Sample test	02/19/00	48	0	0	0	0	0	1kal	5/31/66
EA111113	0	Sample test	09/06/00	248	0	0	0	0	0	2 hbb	8/2/67
EA111114	0	Sample test	09/06/00	248	0	0	0	0	0	3 lsc	5/22/66
EA111115	0	Sample test	09/06/00	248	0	0	0	0	0	4 jmm	4/11/74
EA111117	0	Sample test	09/06/00	248	0	0	0	0	0	5 lmv	9/9/72
EA111116	0	Sample test	09/06/00	98	0	0	0	0	0	1kal	5/31/66
EA111118	0	Sample test	09/06/00	248	0	0	0	0	0	2 hbb	8/2/67
EA111119	0	Sample test	09/06/00	582	0	0	0	0	0	3 lsc	5/22/66
										4 jmm	4/11/74

FIG. 11